

**MULTIDISCIPLINARY SENIOR DESIGN PROJECT  
GE 498**

**COLLEGE OF ENGINEERING  
VALPARAISO UNIVERSITY  
VALPARAISO, INDIANA**

**System Design Requirements  
for  
Biaxial Tensile Tester**

**TEAM TENSILE TESTER  
Date: September 7nd, 2021**

Prepared by:	<b>Team Tensile Tester</b>	
Concurrence:	G. El-Howayek	Initials:
Approved (Customer):	Bethany Luke	Initials:
Approved (Team Leader):	Michael Scarsella	Initials: MS

Total Number of Pages: 5

Revisions (initials in columns indicate approval of revision)						
Letter	Date	Revisor	Advisor	Customer	Team Leader	Reason for Revision
A	9/7	Team	GH	BL	MS	1st revision
B	9/7	Team	GH	BL	MS	Fix requirements
C	9/16	Team	GH	BL	MS	Professor/Customer
D	11/9	Team	GH	BL	MS	Specification Adjust.
E						
F						
G						
H						
I						

**Honor Code Statement**

I have neither given or received, nor have I tolerated other's use of unauthorized aid.

---

Michael Scarsella

---

Garret Brittingham

---

Michael Fermanis

---

David Engstrom

---

Jing Qin

---

Subhanga Ghimire

,

## Goal Statement

The Biaxial Tensile Tester is a device that will test and record the strength of skin, organs and other organic tissues while pulling on the tissue simultaneously in two orthogonal directions.

## Objectives

In order to achieve success, the system shall

1. Be able to pull tissue with the same velocity from all 4 sides simultaneously
2. Be able to measure and record the force and displacement
3. Be able to be set up easily and operated safely
4. Be lightweight and portable

## System Requirements

### 1. TEST REQUIREMENTS

- 1.1 The tissue sample shall be suspended in the air or a fluid during the test to avoid friction on a surface that would distort the measurements.
- 1.2 The system shall apply 0.1 Newtons (in tension) to remove slack from the specimen and threads (if using) before deforming the tissue
- 1.3 The data shall be zeroed after the slack is removed.
- 1.4 The system shall deform tissue in 2 orthogonal directions simultaneously, where the angle between the two directions is within 1 degree of 90.
- 1.5 The tissue shall deform at a constant velocity of 50 +/- 0.1 mm/min in both the x axis and the y axis, which are two orthogonal directions in the horizontal plane.
- 1.6 The system shall maintain an even distribution of forces across each edge of the tissue every 0.1 to 0.5 +/- 0.1 cm each side.
- 1.7 The system shall be able to apply at least 80 Newtons (in tension) of force in both the x direction and the y direction.
- 1.8 The initial size of the specimen shall be between 1 cm by 1 cm and 4 cm by 4 cm.
- 1.9 The maximum positive displacement in both the x and y direction shall be at least 6 cm.

## 2. DATA REQUIREMENTS

- 2.1 The system shall record the actual forces applied within +/- 0.1 Newtons in the x direction and in the y direction.
- 2.2 The system shall record the displacement of the tissue within +/- 0.2 mm in the x direction and in the y direction.
- 2.3 The data shall be exported to a spreadsheet on a computer, SD card, or flash drive.

## 3. SAFETY AND EASE OF OPERATION REQUIREMENTS

- 3.1 [1] The electronics shall have an IPX-1 rating.
- 3.2 The system shall have a start and stop button for collecting data.
- 3.3 The system shall have an emergency stop function.
- 3.4 An untrained student shall be able to learn how to set up the device and perform tests after 30-60 minutes of training.

## 4. PORTABILITY REQUIREMENTS

- 4.1 The system footprint shall be smaller than 3x3ft
- 4.2 The system shall not exceed 6 feet in the vertical direction.
- 4.3 The system shall be powered by a 110V wall outlet.

## Desired Criteria

- 5.1 The velocity may be adjusted to 20 mm/min in the x and y axis.
- 5.2 The velocity may be adjusted to 100 mm/min in the x and y axis.
- 5.3 The velocity may be adjusted to a constant velocity entered by the user (mm/min) in the x and y axis.
- 5.4 The system may include a fluid bath around the specimen at a temperature of 37 degrees C.
- 5.5 The system may record displacement with a camera above, where the frames are synchronized with the force data.

### **Works Cited**

[1] Alam, Noor, 2020, "IPX WaterProof Rating" [Online] Available:

<https://speakersmag.com/ipx-waterproof-rating-ipx4-ipx5-ipx6-ipx7-ipx8-ratings>

[Accessed: 7-September-2021]